



Dosell[®]

Increased safety for patients in need
of long term medication, a tool for
caregivers.

A photograph of an elderly man with white hair, smiling broadly while sitting in a hospital bed. He is wearing a blue hospital gown and holding a white mug. In the background, a female nurse in pink scrubs is smiling and looking down at a tablet. The room is bright and modern, with a white lamp and a white medical device labeled 'Dysell' on the wall. On a bedside table in the foreground, there is a plate of food, a glass of orange juice, and a small bowl of fruit.

A pleasant life for patients,
relatives and healthcare
providers.

Dosell is a patented medicine dispenser for safer medication at home.°

In line with an increasing proportion of older people and a population that consume more and more medicine, the risk of drug injuries escalates. Daily, on average **300** swedes are hospitalized due to incorrect medication. Consequently, individuals suffer unnecessarily and society costs increases.

Dosell eliminates this issue by making sure that correct medicine is taken at the correct time. Additionally, Dosell alerts healthcare personnel and relatives if the medicine sachet is not taken. This is why Dosell reduces the risk of overdosing, taking the wrong combination of medicines or not taking the medication at all.

Dosell is simply the given supplier of safe medication for you as a patient or for your relatives, therewith a valuable tool for you as a healthcare provider.



This is how the medicine dispensing robot works

- 1 Dosell is loaded with the patient's prescribed medication packaged in sachets. One charge contains medicine for two weeks of use.
- 2 When it is time for the patient to take their medicines, Dosell automatically delivers the right sachet at the right time.
- 3 The sachet opens and convey to an open compartment on the front of the appliance.
- 4 When the sachet has fallen down to be taken, a lamp lights up the compartment together with a sound signal reminding the patient that it is time to take the medicine.
- 5 If the patient does not take their medication; the sachet does not leave the compartment, Dosell alerts healthcare personell and relatives who can quickly prevent an injury from happening.



- ✓ Dosell cuts open the medicinal sachet
- ✓ Dosell is able to handle all medical sachets unless it exceeds **160mm** in width
- ✓ Design that is simple and easy to maneuver
- ✓ Dosell enables remote monitoring of medication consumption
- ✓ Dosell is managed through an app available to both iOS and Android
- ✓ Dosell is able to automatically delegate alarm outside office hours to another operator

Dosell increases patient safety
and saves time for healthcare
personnel





Safety for patients and relatives

- Reduced risk of medication errors and injuries
- Ability to live an independent life at home longer
- Increased safety for relatives who can remotely monitor the patient's medication



Increased comfort for staff in health and home care

- Great reduction of medicine handling time
- Saves time and resources for monitoring and controlling the patient's medication
- Opportunity for healthcare personnel to intervene early if lack of medication



Significant savings for the society

- Every reduced home service hour allow savings of SEK 472 *
- Every avoided hospitalization generates a saving of SEK 55,330 **
- 10 % fewer care cases result in a saving of SEK 1.7 million / day - more than half a billion SEK a year

*SKL KPB databasen 2016

**SKL KPP databasen 2015

If you want to know more...

Feel free to contact us or visit our website www.dosell.com



Helena Kolvik

Business unit manager

helena.kolvik@dosell.com

+46 (0)70 540 05 16

Dosell[®]

www.dosell.se